


WIRED 2024



MEG + EEG ANALYSIS & VISUALIZATION



# MNE-Python Tutorial

March 13, 2024  
9:30am – 12:30pm  
Paris, France



Maansi Desai, PhD

Postdoctoral Fellow  
UT Austin



Alexandre Gramfort, PhD

Senior Research Scientist  
Meta Inc.




Apolline Mellot

PhD Candidate  
Univ. Paris-Saclay

# THE AUDITORY PATHWAY



of the human brain



Have you got enough blankets?

hɪv y ɡɑr ɪn ʌf blæŋ kɪts

hɪv y ɡɑr ɪn ʌf blæŋ kɪts



sentence

meaning

word

syllable

phoneme

spectrotemporal  
(frequency decomposition)


acoustic vibration



# INTRACRANIAL RECORDINGS IN THE HOSPITAL




Miller Nature Human Behavior 2019




*Intracranial recordings + cognitive tasks for research*

# HOW DO SIGNALS DIFFER?




Mercier et al. NeuroImage 2022



Chang et al. Neuron 2015


# HIGH GAMMA AS AN INDEX OF NEURAL ACTIVITY





# ADVANTAGES OF INTRACRANIAL RECORDINGS

- Can measure changes at fast timescales
- Phonological information: /b/ vs. /p/
- You know where your signals are coming from
- Very high SNR, requires fewer trials than EEG/MEG



# DISADVANTAGES OF INTRACRANIAL RECORDINGS

- Coverage depends on clinical need, not your research question
- Invasive (must be done in patient populations)
  - Epilepsy/tumor
- Relatively time-limited experiments / less controlled experiments

# CONSIDERATIONS WHEN PREPROCESSING iEEG DATA

- Environment is not as controlled as a typical recording environment (sound proof booth)
- Background sounds
- Where electrodes are located
- How electrodes are plugged into amplifier is important for (re)-referencing
- Sources of electrical noise
- Acquisition sampling rate is important if you want high frequency information







# DIFFERENT TYPES OF iEEG ARTIFACTS

- Motion artifact
- Epileptiform activity
- Chewing/jaw artifact
- VNS artifact (vagus nerve stimulator)

# SEIZURE ACTIVITY (LEFT)



# iEEG TUTORIAL

- We will load iEEG data with BIDS metadata
- Plot power spectrum
- Plot raw data and look for bad channels and artifacts
- Convert to high-gamma analytic amplitude
- Plot evoked data